

# Aperio ImageScope - Manual

## Description

Freely downloadable ImageScope viewing software—experience rapid access to crisp, true-color digital slide images to which you can adjust magnification, pan and zoom, compare different stains, annotate areas of interest, perform image analysis, and more.

Package Notes/Comments – Copy to all **Notes/Comments** sections in SCCM creating package.  
Free software license.

Source file location (share or weblink):

<http://www.leicabiosystems.com/pathology-imaging/epathology/aperio-university/imagescope/>

Copy source files to Software Vault share location:

\\server\share\Sources\Software Vault\Aperio ImageScope

Create a PowerShell script, *Aperio\_ImageScope\_Uninstall.ps1*, to remove Aperio ImageScope

Import *Image-Scope-Rev-v12.1.0.5029.exe* into SCCM Client Applications as a *Script Installer*.

## **Installation program:**

"Image-Scope-Rev-v12.1.0.5029.exe"

## **Uninstall program:**

Powershell.exe –executionpolicy Bypass –file "Aperio\_ImageScope\_Uninstall.ps1"

**Detection Method:** Registry key exists.

Registry hive: HKEY\_LOCAL\_MACHINE

Registry key: SOFTWARE\Wow6432Node\Microsoft\Windows\CurrentVersion\Uninstall\{A5856584-F090-4FD3-BA95-34E6D85546B1}

## **User Experience**

Installation behavior: Install for system

Logon requirement: Whether or not a user is logged on

Installation program visibility: Hidden

Enforce specific behavior: No specific action

**Dependencies** (list redistributables, runtime, etc., available for install with SCCM):

Visual C++ 2010 x64 Redistributable - 10.0.40219

## **Notes**

Installs 2 components:

Aperio ImageScope

GUID: {A5856584-F090-4FD3-BA95-34E6D85546B1} (32-bit)

Registry key or file path (used to check if application is installed):

HKEY\_LOCAL\_MACHINE\SOFTWARE\Wow6432Node\Microsoft\Windows\CurrentVersion\Uninstall\{A5856584-F090-4FD3-BA95-34E6D85546B1}

UninstallString (determines uninstall method to use):

"C:\Program Files (x86)\InstallShield Installation Information\{A5856584-F090-4FD3-BA95-34E6D85546B1}\setup.exe" -runfromtemp -l0x0409 -removeonly

DisplayName: Aperio ImageScope

DisplayVersion: 12.1.0.5029

Publisher: Aperio Technologies, Inc.

#### AperioPrerequisites

GUID: {07138E24-FE31-4AB3-BA41-36A56CDE2815} (32-bit)

Registry key or file path (used to check if application is installed):

HKEY\_LOCAL\_MACHINE\SOFTWARE\Wow6432Node\Microsoft\Windows\CurrentVersion\Uninstall\{07138E24-FE31-4AB3-BA41-36A56CDE2815}

UninstallString (determines uninstall method to use):

MsiExec.exe /I{07138E24-FE31-4AB3-BA41-36A56CDE2815}

DisplayName: AperioPrerequisites

DisplayVersion: 12.00.00.0004

Publisher: Aperio Technologies, Inc.

#### **Aperio\_ImageScope\_Uninstall.ps1**

```
<#
.SYNOPSIS
    SCCM Uninstall program for Aperio ImageScope

.DESCRIPTION
    Remove the following component(s) from window 7 x64 systems:
    - Aperio ImageScope
    - AperioPrerequisites

.NOTES
    FileName: Aperio_ImageScope_Uninstall.ps1
    Author: Jerry Senff
    Updated: 02/28/2019
    Comments: Powershell.exe -executionpolicy bypass -file
"Aperio_ImageScope_Uninstall.ps1"
#>

# Error file
$startLocation = Get-Location
$error=@()
$errorpath = '\ErrorLogs\Aperio_ImageScope_Uninstall.txt'
$errorFileLocation = ${env:SystemDrive} + $errorpath

# Uninstall variables - Aperio ImageScope
$MsiExec = (${env:SystemRoot} + '\System32\msiexec.exe')
$imageScope= 'Aperio ImageScope'
```

```

$ImageScopeUninstall = '"' + ${env:ProgramFiles(x86)} + '\InstallShield Installation
Information\{A5856584-F090-4FD3-BA95-34E6D85546B1}\setup.exe"'
$ImageScopeParams = '-runfromtemp -10x0409 -removeonly'
$ImageScopeRegkey =
'HKLM:SOFTWARE\Wow6432Node\Microsoft\windows\CurrentVersion\Uninstall\{A5856584-F090-
4FD3-BA95-34E6D85546B1}'

# Uninstall variables - AperioPrerequisites
$AperioPrereq = 'AperioPrerequisites'
$AperioPrereqParams = '/x {07138E24-FE31-4AB3-BA41-36A56CDE2815} /qn /norestart'
$AperioPrereqRegkey =
'HKLM:SOFTWARE\Wow6432Node\Microsoft\windows\CurrentVersion\Uninstall\{07138E24-FE31-
4AB3-BA41-36A56CDE2815}'

Write-Host ""
Write-Host
"*****"
Write-Host ""
Write-Host "$ImageScope Uninstaller"
Write-Host ""
Write-Host "Purpose: Remove the following components:"
Write-Host "    - Uninstall $ImageScope"
Write-Host "    - Uninstall $AperioPrereq"
Write-Host ""

# Uninstall Aperio ImageScope
set-location $startLocation
Write-Host
"*****"
Write-Host ""
Write-Host "Checking for $ImageScope installation..."
Write-Host ""

If ((Test-Path -Path $ImageScopeRegkey) -eq 'True') # Check uninstall registry key or
install path
{
    Write-Host "    $ImageScope is installed."
    Write-Host ""
    Write-Host "Uninstalling $ImageScope..."
    Write-Host ""
    Write-Host "Command: $ImageScopeUninstall $ImageScopeParams"
    Write-Host ""
    Start-Process -FilePath $ImageScopeUninstall -ArgumentList $ImageScopeParams -
ErrorVariable +err -Verb Open -wait
    Write-Host ""
    Write-Host "    $ImageScope uninstall complete."
    Write-Host ""
}
Else
{
    Write-Host "    $ImageScope is not installed."
    Write-Host ""
}

start-sleep 5

# Uninstall AperioPrerequisites
set-location $startLocation
Write-Host
"*****"
Write-Host ""
Write-Host "Checking for $AperioPrereq installation..."
Write-Host ""

If ((Test-Path -Path $AperioPrereqRegkey) -eq 'True')
{
    Write-Host "    $AperioPrereq is installed."
    Write-Host ""
    Write-Host "Uninstalling $AperioPrereq..."
    Write-Host ""
    Write-Host "Command: $MsiExec $AperioPrereqParams"
    Write-Host ""
}

```

```
        Start-Process -FilePath $MsiExec -ArgumentList $AperioPrereqParams -ErrorVariable  
+err -Verb Open -wait  
        Write-Host ""  
        Write-Host "    $AperioPrereq uninstall complete."  
        Write-Host ""  
    }  
    Else  
    {  
        Write-Host "    $AperioPrereq not installed."  
        Write-Host ""  
    }  
  
    Start-Sleep 5  
  
    $err | Out-File $errFileLocation  
    Write-Host "*** Log file location = " $errFileLocation " ***"  
  
    Start-Sleep 5
```