

3dMD Patient

Description

3dMDpatient is a software application that allows a surgeon to display 3D surfaces and take measurements from them. This is done by loading one or more surfaces and manipulating them using the various tools supplied in the software.

Package – Copy to all **Notes/Comments** sections in SCCM during package creation.
Free software license.

Source file location (share or weblink):
\\server\share\files\apps\3dMD Patient

Copy source files to Software Vault share location:
\\server\share\Sources\Software Vault\3dMDpatient_413\3dMDpatient_413_x86

Import *VsInstaller.msi* into SCCM Client Applications as a *MSI* Installer

Installation program:

msiexec /i "VsInstaller.msi" /qn /norestart

Uninstall program:

msiexec /x {0D8D003F-B5AB-4D16-92FE-364A0EA3E1DA} /qn /norestart

Detection Method:

MSI Product Code: {0D8D003F-B5AB-4D16-92FE-364A0EA3E1DA}

User Experience

Installation behavior: Install for system
Logon requirement: Whether or not a user is logged on
Installation program visibility: Hidden
Enforce specific behavior: No specific action

Dependencies (list redistributables, runtime, etc., available for install with SCCM):
None

Notes

GUID: {0D8D003F-B5AB-4D16-92FE-364A0EA3E1DA} (32-bit)

Registry key or file path (used to check if application is installed):
HKEY_LOCAL_MACHINE\SOFTWARE\Wow6432Node\Microsoft\Windows\CurrentVersion\Uninstall\{0D8D003F-B5AB-4D16-92FE-364A0EA3E1DA}

UninstallString (determines uninstall method to use):
MsiExec.exe /X{0D8D003F-B5AB-4D16-92FE-364A0EA3E1DA}

DisplayName: 3dMDpatient

DisplayVersion: 4.1.3

Publisher: 3dMD